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scATAC-seq (Single-Cell Assay for Transposase-Accessible Chromatin using sequencing) allows for assessing genome-wide chromatin accessibility at the single-cell level. The technology is based on the use of the hyperactive transposase Tn5 which cuts exposed open chromatin and simultaneously ligates adapters for subsequent amplification and sequencing. Coupled to BioRad ddSEQ Single-Cell Isolator, this technology allows genome-wide open chromatin sequencing of 500 - 5000 droplet isolated single nuclei.

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In order to maximize results, we provide strict guidelines for sample preparation (cell collection).

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We request 200.000 - 1.000.000 cells / sample of interest to be processed. Cell viability should be higher than 95% when collected. Two additional test samples (200.000 – 1.000.000 cells / test sample) are required to perform the nuclei isolation optimization.

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scATAC-seq workflow

scATAC-seq 解析内容

Cells will be prepared by the customer according to Diagenode's sample preparation guidelines. Diagenode recommends to perform the experiment at least on two biological replicates.

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OPTIMIZATION ON TEST SAMPLES

(2 test samples requested per cell line/tissue type)

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1. Sample QC

- Cell counting
- Cell viability assay

2. Nuclei isolation

- Nuclei isolation: testing 2 lysis conditions
- QC: Nuclei counting and microscopic evaluation

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1. Sample QC

- Cell counting
- Cell viability assay

2. Sample tagmentation

- Nuclei isolation
- QC: Nuclei counting and microscopic evaluation
- Tn5 digestion

3. Droplet encapsulation and library preparation

- Isolate single nuclei by encapsulation
- Library preparation: single nuclei barcoding, sample indexing and amplification
- Library purification
- QC on the scATAC-seq library: DNA concentration and analysis of library profile
- Libraries pooling

4. Deep sequencing

- Samples are sequenced on an Illumina® platform.
- Sequencing depth will be adjusted to project dependent criteria (25-50K raw reads/cell on average).

5. Bioinformatic analysis

- Quality check.
- Debarcoding.
- Alignment to reference genome.
- Bead filtration
- Bead deconvolution
- Identification of enriched regions (peak calling).
- Generation of cells-by-peaks count matrix.

Provided files:

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- Reports with sequencing statistics and cell encapsulation efficiency.
- Raw data in FASTQ format.
- FastQC reports.
- Alignment files in BAM format.
- Peak files in BED format.
- Count matrix in csv and sparse mtx format.

GUIDELINES

Additional analysis on request: 追加解析オプション

- **Cell classification:** Identification of cells subpopulations by unsupervised clustering (using Seurat/Signac).
- **Differential accessibility analysis:** Identification and annotation (human, mouse, rat, drosophila) of differentially accessible regions between cell subpopulations.
- **Gene ontology terms analysis:** Enrichment analysis on genes associated with differentially accessible regions between cell subpopulations. Gene Ontology terms that are overrepresented in differentially accessible regions may indicate the underlying biological processes involved.
- **Pathway analysis:** Identify biochemical pathways in which genes associated with differentially accessible regions between cell subpopulations may be overrepresented.
- **Visualization of specific genomic regions:** Visualization of results (i.e. sequencing data, peaks) at specific genomic regions (e.g. genes, promoters) for each cell subpopulations in publication-ready images (human, mouse, rat).

Additional information

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For sample preparation and sample shipment it is mandatory to follow Diagenode's guidelines. If customer samples do not meet Diagenode's quality requirements, any additional QC of new samples will be charged to the customer. Any delay in sample shipment to Diagenode's facilities might result in delaying customer's project.

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Generated files will be available for download during 1 month and stored for an additional period of 3 months on Diagenode's servers. Additional long-term storage of data is available upon request. This offer includes a one hour call to walk you through the results if needed.

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Original samples are stored at Diagenode during 12 months after project completion, but will be discarded once this time is exceeded. Return shipment of samples is available upon request.

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Any additional service which is beyond the current project scope will be charged to the customer.

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