




MeDIP-qPCR or MeDIP-Seq?

MeDIP-qPCR	MeDIP-Seq
single-locus data requires design of loci-specific primers	genome-wide data suitable for discovery analysis
low-cost fast	high sequencing cost longer workflow (library prep, sequencing, data analysis)



Workflow: Methylated-DNA ImmunoPrecipitation (MeDIP)





Cell lysis & gDNA extraction



High
integrity gDNA

- **Cultured cells**

Recommended kit: [XL GenDNA Extraction Module](#) (Diagenode)


For large amounts of cultured cells (~1-1.5 million)

- **Other template (e.g. tissue, blood)**

Flexibility to use other kits

Ensure good gDNA quality and purity


- **Starting amount for MeDIP: 10ng - 1 µg of sheared gDNA**




Whole Blood



Fresh Tissue



Cell



FFPE Tissue



gDNA shearing

- Shear gDNA to consistent fragment mean size
 - 300-400 bp for qPCR
 - 200 bp for NGS
- Use a good sonicator
 - Recommended: [Bioruptor® Pico](#) (Diagenode)
- Shearing TE buffer
- Sample concentration at 2-100 ng/μl



Precise
fragmentation



gDNA shearing: Parameters for Bioruptor® Pico



Guide for successful gDNA shearing using the Bioruptor® Pico

Sample type	Tube holder	Tube	Volume of sample	Shearing protocol
Low sample amounts (from 10 ng)	for 0.2ml tubes B01201144	0.2ml Bioruptor® Microtubes C30010020	50 µl	13 cycles [30sec ON / 30sec OFF]
Standard sample amounts	for 0.65ml tubes B01201143	0.65ml Bioruptor® Microtubes C30010011	100 µl	13 cycles [30sec ON / 30sec OFF]

<https://www.diagenode.com/files/protocols/protocol-dna-shearing-bioruptor-pico.pdf>


Features & Benefits

- Best suited for NGS sample preparation
- Superior sample yields
- Consistent fragment size
- High reproducibility




gDNA shearing: Examples of DNA profiles

BEFORE



AFTER



Performing IP



Keep on ice to
maintain ssDNA

1. Prepare Beads (washing and resuspension)




2. Prepare IP reaction tube

Sheared gDNA
Spike-in controls
Buffers

3. heat-denaturation to ssDNA




4. Set aside INPUT samples



5. Add 5mC specific antibody

Overnight incubation on rotating wheel



6. Remove unmethylated DNA

Wash bead/antibody/DNA complexes



Performing IP: INPUT sample

- Fraction of sheared DNA is kept aside as INPUT (~10%)
Processed in parallel with IP-samples (without antibody)
Recommended: 1 input for each DNA sample
- Key reference for MeDIP-qPCR and MeDIP-seq analysis
- MeDIP-qPCR: used to calculate the IP recovery (% of input)
- MeDIP-seq: mandatory for bio-informatics analysis
Normalization for mappability of a region, to avoid duplication bias etc.



me-DNA IP
samples +
INPUT controls





Performing IP: 5-mC antibody

Recommended antibody:

5-methylcytosine (5-mC) Antibody - clone 33D3 (Diagenode)



Highly specific
antibody





Performing IP: Spike-in controls


- **External DNA spike-in controls** & primer pairs
 - Synthetic sequence without homology to any model species
 - No interference with DNA sample of interest
 - Add directly to each sheared gDNA sample before IP
 - Not sequenced but use as control during qPCR analysis
 - **Methylated (positive) spike-in:** confirm overall efficiency of IP workflow
 - **Unmethylated (negative) spike-in:** measure non-specific IP background
- Recommended kit: [DNA methylation control package V2](#) (Diagenode)



External spike-
in controls



DNA isolation



qPCR analysis: QC check

IP QC check using qPCR on the spike-in controls

1. % recovery

$$= 2^{[Ct(INPUT) - 3.32 - Ct(IP)]} * 100$$



$= \log 10 / \log 2$, where 10 is dilution factor IP/INPUT

2. Enrichment


$$= \frac{\% \text{ recovery meDNA}}{\% \text{ recovery unDNA}}$$

3. Specificity

$$= 1 - \frac{1}{\text{Enrichment}}$$

Success criteria

- Methylated spike-in recovery > 20%
- Unmethylated spike-in recovery < 1%
- Spike-in IP specificity > 99%





qPCR analysis: QC check & analysis




qPCR quality controls

- Primer pairs for **Positive/Negative control regions** within the DNA sample of interest


QC check during qPCR & NGS analyses

- Available [primer pairs](#) for human/mouse/rat
 - Positive TSH2B locus:** methylated region
 - Negative GAPDH locus:** unmethylated region
- Or design your own qPCR primers






Summary – Tips for Good MeDIP-qPCR Assay




Workflow: Methylated DNA ImmunoPrecipitation (MeDIP)





Library preparation: Adaptor ligation before IP

- Critical step for successful sequencing
- Adaptor ligation on dsDNA before denaturation and IP
- Recommended kit: [iDeal DNA IP Library Preparation Kit](#)
 - Fully optimized & validated for MeDIP-seq
 - Low input (from 10ng DNA)
 - High library yield



MeDIP-seq: Performing IP



Keep on ice to
maintain ssDNA

1. Prepare Beads (washing and resuspension)




2. Prepare IP reaction tube

Adaptor-ligated DNA
Spike-in controls
Buffers

3. Heat-denaturation to ssDNA




4. Set aside INPUT samples



5. Add 5mC specific antibody

Overnight incubation on rotating wheel




6. Remove unmethylated DNA

Wash bead/antibody/DNA complexes



MeDIP-seq: DNA isolation





Library preparation: Optimal PCR cycles

Minimization of PCR cycles

- Limits the risk of introducing bias during PCR
- Minimizes # of PCR duplicates




Minimize #
PCR cycles

Determine optimal # of PCR cycles for each sample and input on a small aliquot

- **Optimal cycle number** = Ct (*rounded up*) + 1


e.g. $Ct = 6.82 \rightarrow 8$ amplification cycles





Library preparation: Amplification

- Recommended kit: [iDeal Unique Dual Indexes for Illumina](#)
 - Library multiplexing
 - Unique i5 and i7 barcodes
 - Filtering of index-hopping events (especially on NovaSeq)





Library preparation: Pooling & quality controls

- **Determine library size**


- Bioanalyzer or Fragment Analyzer

- **Quantify**

- Qubit
- Convert ng/ μ l to nM by average library size

- **Dilute and pool normalized libraries**

- Same size for best clustering





MeDIP-seq: Sequencing Settings

Read length

- Paired-end 50 bp

Sequencing depth

- Mainly set by genome size, samples and flow cell type
- 50 M for genome-wide DNA methylation analysis (human & mouse)
- Use same depth for input

Replicates

- Biological duplicates are minimum for most bioinformatic analyses
- Increased replicate number will improve sensitivity of the downstream analysis

Input sequencing

- One input per sample is gold standard
- Pooling inputs from replicates can often be considered






MeDIP-seq: Analysis

Bioinformatician	R	https://bioconductor.org/packages/release/bioc/html/MEDIPS.html
	free-ware and online tool kits	www.bioinformatics.babraham.ac.uk/projects/seqmonk/ biit.cs.ut.ee/gprofiler/gost
Wet-lab expert with free time	standard bioinformatic services	https://www.diagenode.com/en/categories/Services Comprehensive DNA methylation and bioinfo services
	free-ware and online tool kits	www.bioinformatics.babraham.ac.uk/projects/seqmonk/ Initial & advanced data analysis, genome browser, graphical presentation of data https://biit.cs.ut.ee/gprofiler/gost Functional profiling tool
No expertise or no free time	advanced bioinformatic services	https://www.diagenode.com/en/categories/Services Comprehensive DNA methylation and bioinfo services




MeDIP-seq: Analysis





MeDIP-seq: Example of results

Sequencing Statistics



Uniquely mapped

= reads after mapping

+ duplicates removal

+ multi-mappers removal

= reads for methylation analysis


Sample: Human whole blood genomic DNA

Starting amount: 1 µg; 100 ng; 50 ng; 25 ng; 10 ng

Sequencing: PE 50 bp on Illumina Novaseq

Bioinformatics analysis: MeDIPS


MeDIP-seq: Enrichment visualization



Enrichment of **me-DNA** compared to **INPUT** sample




MeDIP-seq: Coverage profiles






MeDIP-seq: Peak calling






Summary – Tips for Good MeDIP-seq Assay





Summary – Diagenode Support




Hardware:	<u>Bioruptor</u>			
Reagents:	<u>XL GenDNA Extraction Module</u>	<u>5-mC antibody Spike-in controls</u>	<u>IPure kit</u>	<u>Primer Pairs</u>
Kits:	<u>MagMeDIP qPCR Kit</u>			
	<u>MagMeDIP-seq Package V2</u>			
Services:	<u>Epigenomic Profiling Services</u>			
	<u>Data Analysis Service</u>			

Applications

MeDIP-qPCR on FFPE samples

Mol Neurobiol
DOI 10.1007/s12035-016-0345-x

Evaluating the Feasibility of DNA Methylation Analysis Using Long-Term Archived Brain Formalin-Fixed Paraffin-Embedded Samples

Stine T. Bak¹ · Nicklas H. Staunstrup^{1,2,3,4} · Anna Starnawska^{1,3,4} · Tina F. Daugaard¹ · Jens R. Nyegaard^{5,6} · Mette Nyegaard^{1,3} · Anders Borglum^{1,3,4,7} · Ole Mors^{3,4,7} · Karl-Anton Dorph-Petersen^{2,6,8} · Anders L. Nielsen^{1,3} 

Received: 3 August 2016 / Accepted: 5 December 2016
© Springer Science+Business Media New York 2016

MeDIP-seq on cfDNA from plasma...

LETTER

<https://doi.org/10.1038/s41586-018-0703-0>

Sensitive tumour detection and classification using plasma cell-free DNA methylomes

Shu Yi Shen^{1,12}, Rajat Singhania^{1,12}, Gordon Fehringer^{2,12}, Ankur Chakravarty^{1,12}, Michael H. A. Roehri^{1,12,14}, Dianne Chadwick¹, Philip C. Zuzarte¹, Ayelet Borgida², Ting Ting Wang^{1,14}, Tian Tian Li¹, Olena Kis¹, Zhen Zhao¹, Anna Spreafico¹, Tiago da Silva Medina¹, Yadan Wang¹, David Roulois^{1,6}, Ilias Eitayebi^{1,4}, Zhuo Chen¹, Signy Chow¹, Tracy Murphy¹, Andrea Arruda¹, Grainne M. O'Kane¹, Jessica Liu¹, Mark Mansour¹, John D. McPherson⁷, Catherine O'Brien¹, Natasha Leigh¹, Philippe L. Bedard¹, Neil Flesher¹, Geoffrey Liu^{1,14}, Mark D. Minden¹, Steven Gallinger^{1,10}, Anna Goldenberg¹, Trevor J. Fughi^{1,4}, Michael M. Hoffman^{1,13}, Scott V. Bratman^{1,4}, Rayjean J. Hung^{2,18} & Daniel D. De Carvalho^{1,4*}

MeDIP-qPCR on cfDNA

IMCM
Winter 2017, Vol 6, No 1

MeDIP Real-Time qPCR has the Potential for Noninvasive Prenatal Screening of Fetal Trisomy 21

Mohammad Kazemi^{1,2,3}, Mansoor Salehi^{1,2,3*}, Majid Kheirollahi¹

1. Department of Genetics and Molecular Biology, School of Medicine, Isfahan University of Medical Sciences, Isfahan, Iran.
2. Medical Genetic Center of Genome, Isfahan, Iran.
3. Pediatric Inherited Diseases Research Center, Research Institute for Primordial Prevention of Noncommunicable Disease, Isfahan University of Medical Sciences, Isfahan, Iran.

MeDIP-seq on plants



The methylation pattern of DNA and complex correlations with gene expressions during *TuMV* infection in Chinese cabbage

J. YU¹, L.-W. GAO¹, Y. YANG¹, C. LIU¹, R.-J. ZHANG¹, F.-F. SUN², L.-X. SONG³, D. XIAO¹, T.-K. LIU¹, X.-L. HOU¹, and C.-W. ZHANG^{1*}

State Key Laboratory of Crop Genetics and Germplasm Enhancement, College of Horticulture, Nanjing Agricultural University, Nanjing, 210095, P.R. China¹
Nanjing Vegetable Science Research Institute, Nanjing, 210095, P.R. China²
Jiangsu Key Laboratory for Horticultural Crop Genetic Improvement, Jiangsu Academy of Agricultural Sciences, Nanjing, Jiangsu, 210014, P.R. China³

MeDIP-seq on oysters



RESEARCH ARTICLE

Dynamics of DNA methylomes underlie oyster development

Guillaume Riviere^{1,2*}, Yan He^{2*}, Samuele Tecchio^{1,2*}, Elizabeth Crowell^{1,2}, Michael Gras^{1,2}, Pascal Sourdis^{1,2}, Ximing Guo¹, Pascal Favre^{1,2}

¹ Normandy University, Caen, France, ² Université de Caen Normandie, UMR BOREA MNHN, UPMC, UCBN, CNRS-7208, IRD-207, Caen, France, ³ Ministry of Education Key Laboratory of Marine Genetics and Breeding, College of Marine Life Sciences, Ocean University of China, Qingdao, Shandong, China, ⁴ Haskin Shellfish Research Laboratory, Department of Marine and Coastal Sciences, Rutgers University, Port Norris, N.J., United States of America

* These authors contributed equally to this work.
* guillaume.riviere@unicaen.fr

... and urine



BRIEF COMMUNICATION

<https://doi.org/10.1038/s41591-020-0933-1>



Detection of renal cell carcinoma using plasma and urine cell-free DNA methylomes

Pier Vitale Nuzzo^{1,2,10}, Jacob E. Berchuck^{1,2,10}, Keegan Korthauer^{1,10}, Sander Spisak^{2,10}, Amin H. Nassar^{1,2}, Sarah Abou Alaiw^{1,2}, Ankur Chakravarty¹, Shu Yi Shen¹, Ziad Bakouny¹, Francesco Boccardo^{1,7}, John Steinharter¹, Gabrielle Bouchard¹, Catherine R. Curran¹, Wenting Pan¹, Sylvan C. Bacal^{1,2}, Ji-Heui Seo^{1,2}, Gwo-Shu Mary Lee^{1,2}, M. Dror Michaelson^{1,9}, Steven L. Chang¹⁰, Sushrut S. Waikar¹⁰, Guru Sonpavde¹, Rafael A. Irizarry^{10,11}, Mark Pomerantz^{1,2}, Daniel D. De Carvalho^{1,10,17}, Toni K. Choueiri^{1,10,17} and Matthew L. Freedman^{1,2,10,17,22}