

Sample Submission Guidelines for RNA-seq

Contact: services@diagenode.com

Website: <https://www.diagenode.com/en/categories/rna-seq-category>

1. Preparation of the samples

The quality of the RNA to be used for RNA-seq is very important and depending on your application we will send you the adequate RNA sample preparation guidelines.

2. Shipment Instructions

- Prepare the samples in 1.5 mL or 0.5mL tubes and label them carefully.
- Fill the **Sample_RNA_ID_sheet.xlsx** and the **Biosafety_form.doc** file with all relevant information and send them by email to **services@diagenode.com**.
- Always send the samples on dry ice and make sure you supply sufficient dry ice for the whole transport and unexpected delays.
- Print the **Sample_RNA_ID_sheet.xlsx** and the **Biosafety_form.doc** and include them in your package when shipping the samples.
- Attach a Dry Ice Sticker with the net weight (see below) on the outside of the parcel. Please send the parcel at the beginning of the week so it gets to us before the weekend.



- Please contact your courier (DHL, TNT, FedEx, UPS, etc) to arrange the shipment of the samples to the following address:

Catherine Creppe
Diagenode
RNA-seq service
Liège Science Park
Rue du Bois Saint-Jean, 3
4102 Seraing (Ougrée)
Belgium

- Communicate the tracking number to **services@diagenode.com**.